






VB*EMS ELASTOMERIC RESPIRATOR SUPPORTED COMBINATIONS			REUSABLE if exposed	After exposure	Comment
	May use w/ HALF facepiece	May use w/ FULL facepiece			
3M 7093 P100 cartridges (magenta) 	Yes 1st choice	Yes 1st choice	Yes	Wipe outside w/ non-bleach disinfectant, keep inside dry	"7093"
3M 7093C P100/HF cartridges (magenta/olive) 	Yes 2nd choice	Yes 2nd choice	Yes	Wipe outside w/ non-bleach disinfectant, keep inside dry	Contingency substitute for 7093
3M 2091 P100 filters (magenta) 	Yes 3rd choice	Yes 3rd choice	NO	Discard	"Pink pancakes"
3M 60923 OV/AG/P100 cartridges (magenta/yellow) 	NO	ONLY COMBO FOR CIVIL UNREST	Yes	Wipe outside w/ non-bleach disinfectant, keep inside dry	

NOTE: Hard cartridges are preferred because only hard cartridges can be wiped off and reused after exposure, and/or used in the rain.